

Permanent Makeup Technology

Internationally Accredited




Product Catalogue

2010





10% Student discount on all purchases for all Permanent Makeup Technology students only.

Rotary Pen Accessories



Black Giant Sun		One of the original manufactured rotary machines. A Wonderful durable machine. Kit comes with rotary machine, ac adaptor, battery pack, caps, needles, sleeves & tubes.	R 1500.00
Giant Sun Tubes		Holds needle tight in the machine. Should be changed after every use.	R 200.00 (10 pack)
Sponges		Inserted on the back end of the needle to prevent any pigment to seep into the motor. Replace every use.	R 50.00 per pack
Giant Sun Covers or Caps LARGE		10 pack 5 to 7 Point LARGE Needle Covers BLACK GIANT SUN	R 50.00
Giant Sun Covers or Caps SMALL		10pack 1 to 3 Prong Small Needle Covers BLACK GIANT SUN	R 50.00
Giant Sun Needles 3-prong		10 Pack 3-Prong BLACK GIANT SUN	R 160.00
Giant Sun Needles 5-prong		10 Pack 5-Prong BLACK GIANT SUN	R 190.00
Giant Sun Needles 7-prong		5Pack 7-Prong BLACK GIANT SUN	R 110.00

Non Iron Oxide Pigments

<p style="text-align: center;">Small Colour Collection</p>		<p>16 Colour collections in 1/6 ounce. Free Colour Mixing guide and Artist colour Wheel. Non Iron Oxide, no colour changes. Get 16 colours for the price of 15.</p>	<p style="text-align: center;">R 1400.00 Full set Save R120.00</p>
<p style="text-align: center;">Medium Colour Collection</p>		<p>16 Colour Collection in 1/3 ounce. Free Colour Mixing Guide and Artist Colour wheel. NON Iron Oxide, no colour changes. Get 16 colours for the price of 15.</p>	<p style="text-align: center;">R 2350.00 Full set Save R210.00</p>
<p style="text-align: center;">Small size Pigments</p>		<p>1/6 ounce NON IRON OXIDE pigment 16 colours to choose from. Refer to colour chart for colour choices.</p>	<p style="text-align: center;">R 95.00 each</p>
<p style="text-align: center;">Medium Size Pigments</p>		<p>1/3 ounce NON IRON OXIDE pigments. 16 colours to choose from. Refer to colour chart for colour choices</p>	<p style="text-align: center;">R 160.00 each</p>

Useful Tips on Non Iron-Oxide Pigments

CAUTION

With anyone's pigments you want to make sure you clean the skin very well before doing a procedure, skin always has bacteria on it. If you open the skin there is a chance of bacterial infection if not cleaned well.


Definition of Iron Oxide is: **Oxygenated Iron.**

Using terms like "inorganic" is other company's ways to confuse the issue to promote their product.

Pigment Removal Solution

Lighten or remove implanted pigment. This product is used to remove or correct permanent makeup or tattoos. Can be used with manual hand method or the machine method

This product is only sold to professionals.

Pigment Removal Solution		Permanent Makeup & Tattoo Removal (10mm) Lighten or remove implanted pigment. Allow 30 days for healing to see results. May need more than one application to remove all colours. Do not get into eyes. Do not use in the eye area.	R 800.00
--------------------------	---	--	-----------------

Permanent Makeup Removal or Correction

Great for taking out old "iron oxide" pigments that have turned a funny colour.

Fix those strange "iron oxide" coloured brows the right way.

Implanting a corrective pigment over an old strange iron oxide pigment will not be a permanent fix....only temporary

1. Take out first. (Applied like implanting)
2. Let heal for 3 months
3. More than one application may be needed
4. Implant a non iron oxide pigment after skin completely heals.

DO NOT OVER WORK THE SKIN

Tattoo Removal


May take more than one session to achieve the desired result.

Always let the skin heal between sessions.

Cleanliness is very important during the healing process; educate your clients well on after care.

**THE PRODUCT WILL WORK FOR YOU. DO NOT DIG WITH YOUR NEEDLES,
YOU MAY CAUSE SCARING.**

Topical Anaesthetic

Blue Gel		Blue Gel is great for open skin procedures contains 4% lidocaine and 2% tetracaine and a vascular constrictor to reduce swelling, bleeding	R640.00
Zone 1		Microcaine 5% Lidocaine provides excellent pain relief for approximately 1 hour when used as directed. Gentle pH of 7.5 (same as human tears.)	R440.00

Permanent Makeup Accessories

Surgical Marking Pens		Don't even worry about losing your design again. Trace around your design.	R 50.00 each
Pigment Cups		Clear Pigment cups 2 Sizes (SMALL or LARGE) Never reuse. Clear helps you see the colour pigment better. Great for Topical, mixing lash & brow tint.	R 65.00 Large 100's R 65.00 Small 100's
Pigment Cup Holder		Stainless Steel Pigment cup holder. Holds 4 small cups and 2 large cups	R 100.00 each
Sharps® container		1-Quart containers feature a translucent lid for easy ID of contents and Dual opening accommodates a variety of sharps	R 100.00 each
After care ointment		25 count. Tube with antibiotic effect.	R 450.00 Per box R18.00 each
Micro Brushers		Micro-applicators Perfect to apply Topical Anaesthetic to very small areas. Use one Micro-applicator per procedure. Dispose after completion of procedure	R 100.00 (50) R190.00 (100)

Handles / Holders






Handles in a variety of styles for a comfortable and EZ to fit your hand.

The Original Yellow & Black Handle		The "Original" Yellow & Black Autoclavable handle holds the slant or round needle configuration. Needle sizes #8 to #78.	R 160.00
Use & Toss 10 Point Very Popular		One time use tool. No assembly required. Needle diameter: 0.35mm	R 60.00 each

Manual Needles



Our precision crafted surgical steel needles produces the most gentle and comfortable touch when applying permanent makeup. Less trauma to the skin will speed up the healing process. Needles are sterilized and should only be used one time only. Never reuse needles. Always discard used needles in a sharps container.

8 Point Needle		8-point needle (0.35) fits into manual handle. Great for touch-ups, eyelash enhancement, and thin brows. Black & Yellow handle or metal handle.	R 50.00
10 Point Needle		10-point needle (0.35mm) fits into the manual handle. Great for touch-ups, eyelash enhancement, and thin brows. Black & Yellow handle or metal handle.	R 50.00
12 Point Needle		12-point needle (0.35mm) fits into the manual handle. Great for touch-ups, eyelash enhancement, eyeliner and brows. Black & Yellow handle or metal handle.	R 50.00
12 Point FINE Needle		12-point needle (0.25mm) FINE Create a perfect hair stroke Black & Yellow handle or metal handle	R 50.00
14 Point Needle		14-point needle (0.35 mm) Good for eyebrows, eyeliner, lip liner Black & Yellow handle or metal handle	R 50.00

